



Domus 4


Alternating Pressure Redistribution System with Low Air Loss

- Suitable for patients with medium to high risk of pressure ulcers
- Ventilated low air loss for patients' utmost comfort
- Cell-on-cell bladder to avoid bottoming out
- Visual and audible alarm enhances patients safety and security




Order Product Code: **9C077B66**
Order Online: **dearjane.com.au**

APEX
Wellell Group




Domus 4

Alternating Pressure Redistribution System with Low Air Loss



Ventilated Low Air Loss
Wicks moisture away to keep the patient comfortable



Cell-on-Cell Bladder
The bottom cells acts as basic pressure support to avoid bottoming out in case of power failure



Zipper-all-Around Cover
Reduce ingress of unwanted fluid and simply one step to remove for cleaning



CPR Valve
Quickly deflate the mattress during emergencies



Nylon/PU Cover
Vapor permeable and water resistant with basic flame retardant level



Biocompatibility



Vapour Permeable



Water Resistance



Flame Resistance



Quick Connector
Once with cap on, instantly switch to static air mattress for transportation need



Visual and Audible Alarm

- Low Pressure
- Power Failure
- Tech Support