



Etac Rufus bath chair

Designed to assist users and carers with every day routines.

 **etac**[®]

Creating Possibilities

Etac Rufus for the bath tub

Rufus bath products make sitting in the bath both easier and safer. To facilitate transfer to and from the bathtub they can be combined with Fresh or Rufus Plus bath board.



Rufus bath stool

The Rufus bath stool allows the user to sit comfortably 15 cm 6" below the edge of the bathtub.



Rufus bath chair

The Rufus bath chair allows the user to sit comfortably and securely 19 cm below the edge of the bathtub. It has a back support to ensure the best comfort and safety. The bath chair is available in two widths

Description

Rufus Bath Chair 70 cm	3.2 kg 7.1 lbs	81506011
Rufus Bath Chair 74 cm	3.2 kg 7.1 lbs	81506012
Rufus bath Stool 70 cm	2.6 kg 5.7 lbs	81504011

Weight

Item no.

Dimensions

Seat area:
41x40 cm 16¼ x 15¾"
External width:
70 x 30.4 cm 27¼ x 12"

Material

Frame: Powder-painted steel tubing
Seat: Polypropylene with embedded aluminium tubing



max 85°C / max 185°F
for maximum of 3 minutes

max 130 kg / max 287 lbs



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.




red dot





 **etac**[®]
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”

 = Standard Stock Item