

iGo

Portable oxygen concentrator





DeVilbiss Healthcare introduces the iGo portable oxygen concentrator.

Small and lightweight, yet capable of delivering up to 3lpm of constant flow oxygen meeting the oxygen requirements for most users.

An ideal choice for travelling, the iGo can be operated from a standard 12 V car charger socket or from its lithium ion battery which will give the user up to 5.4 hours of usage (when set at 1 lpm PulseDose mode). The iGo comes complete with a battery, AC power adapter, mains lead and DC power adapter. The iGo can be transported using a wheeled carry case for added portability. The iGo includes the DeVilbiss Oxygen Sensing Device and service alerts to give piece of mind that the iGo is delivering exactly what is required.

A simple, intuitive control panel makes the iGo easy to use and understand. The iGo is robustly designed to deal with the rigours of today's travel. Other

than cleaning the gross particle filter, preventative maintenance is reduced as the iGo does not require manufacturer's servicing for the first three years.

Main Features

- Up to 3 lpm constant flow
- Up to 6 lpm pulsed flow (equivalent delivery)
- Long battery life
- Oxygen Sensing Device and audible and visual alarms
- Robust design
- Wheeled carry case available
- 12 V car adapter

Typical Battery Operation Times

Flow Rate	Cont. Flow	PulseDose (20 BPM)
1.0	4.0 hrs	5.4 hrs
2.0	2.4 hrs	4.7 hrs
3.0	1.6 hrs	4.0 hrs
4.0	–	3.5 hrs
5.0	–	3.2 hrs
6.0	–	3.0 hrs

SPECIFICATIONS

Dimensions (H x W x D)

38,1 x 27,9 x 20,3 cm

Weight with/without battery

8,6 kg / 7,0 kg

Maximum Altitude

0 - 4000 metres

Nominal Sound Level

40 dB(A) at 3.0 lpm in PulseDose delivery mode

Oxygen Purity

91% ± 3%

Maximum Constant Flow

3 lpm

Maximum PulseDose Flow

6 lpm (equivalent delivery)

Calculated Bolus Size

14 cc per l/min

Alarm Functions


Low oxygen, No flow,
No breath w/failsafe,
Low battery, No power,
Unit malfunction and
High temperature



iGo and wheeled carry case



Wheeled carry case showing extendable handle



iGo Control Panel



iGo Lithium ion Battery

DJMed Pty Ltd
Suite 10/20 Barcoo St
Roseville, NSW 2069
0421 115 992
www.dearjane.com.au

DJMed Pty Ltd
10/20 Barcoo Street,
Roseville, NSW, 2069

web DearJane.com.au
email info@dearjane.com.au
fax (02) 4744 2420
phone 0421 115 992

nsw (02) 8317 5097
qld (07) 3067 7073
vic (03) 8657 7678
sa (08) 7111 2981
wa (08) 6388 8026